


Product datasheet for **TA324618S**

ZNF71 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	This ZNF71 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 74-101 amino acids from the N-terminal region of human ZNF71.
Formulation:	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Purification:	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	54498 Da
Gene Name:	zinc finger protein 71
Database Link:	NP_067039 Entrez Gene 58491 Human Q9NQZ8
Synonyms:	EZFIT
Protein Families:	Druggable Genome, Transcription Factors


[View online »](#)

Product images:


ZNF71 Antibody (N-term) (Cat. #[TA324618])
western blot analysis in human placenta tissue
lysates (35ug/lane). This demonstrates the ZNF71
antibody detected the ZNF71 protein (arrow).