

Product datasheet for **TA324643S**


ZNF307 (ZKSCAN4) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	This ZKSCAN4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 446-473 amino acids from the C-terminal region of human ZKSCAN4.
Formulation:	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Purification:	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	61579 Da
Gene Name:	zinc finger with KRAB and SCAN domains 4
Database Link:	NP_061983 Entrez Gene 387032 Human Q969J2
Synonyms:	P1P373C6; ZNF307; ZNF427; ZSCAN36
Protein Families:	Transcription Factors



[View online »](#)

Product images:

ZKSCAN4 Antibody (C-term) (Cat. #TA324643])
western blot analysis in A549 cell line lysates
(35ug/lane). This demonstrates the ZKSCAN4
antibody detected the ZKSCAN4 protein (arrow).