

Product datasheet for **TA324652S**


TFDP3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	This TFDP3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 82-108 amino acids from the N-terminal region of human TFDP3.
Formulation:	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Purification:	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	44967 Da
Gene Name:	transcription factor Dp family member 3
Database Link:	NP_057605 Entrez Gene 51270 Human Q5H9I0
Synonyms:	CT30; DP4; HCA661



[View online »](#)

Product images:

TFDP3 Antibody (N-term) (Cat. #[TA324652])
western blot analysis in WiDr cell line lysates
(35ug/lane). This demonstrates the TFDP3
antibody detected the TFDP3 protein (arrow).