

Product datasheet for TA324703

Tau (MAPT) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IF, IHC, WB

Recommended Dilution: WB: 1:1000, IHC: 1:50~100, IF: 1:10~50 **Reactivity:** Human (Predicted: Mouse, Rat, Bovine)

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: This TAU antibody is generated from rabbits immunized with a KLH conjugated synthetic

peptide between 677-707 amino acids from the C-terminal region of human TAU.

Formulation: Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.

Concentration: lot specific

Purification: This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by

dialysis against PBS.

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 78928 Da

Gene Name: microtubule associated protein tau

Database Link: NP 058525

Entrez Gene 17762 MouseEntrez Gene 29477 RatEntrez Gene 4137 Human

P10636

Synonyms: DDPAC; FTDP-17; MAPTL; MSTD; MTBT1; MTBT2; PPND; PPP1R103; TAU

Protein Families: Druggable Genome

Protein Pathways: Alzheimer's disease, MAPK signaling pathway




OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn


Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Product images:




Western blot analysis of TAU (arrow) using rabbit polyclonal TAU Antibody



Formalin-fixed and paraffin-embedded human cancer tissue reacted with the primary antibody, which was peroxidase-conjugated to the secondary antibody, followed by AEC staining. This data demonstrates the use of this antibody for immunohistochemistry; clinical relevance has not been evaluated. BC = breast carcinoma; HC = hepatocarcinoma.





Confocal immunofluorescent analysis of TAU Antibody (C-term) (Cat#TA324703) with MCF-7 cell followed by Alexa Fluor 488-conjugated goat anti-rabbit IgG (green).DAPI was used to stain the cell nuclear (blue).