

Product datasheet for **TA324869**


BAIAP2 Mouse Monoclonal Antibody [Clone ID: 1037CT15.1.3]

Product data:

Product Type:	Primary Antibodies
Clone Name:	1037CT15.1.3
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human, Mouse
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Purified His-tagged BAIAP2 protein was used to produced this monoclonal antibody.
Formulation:	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Purification:	This antibody is purified through a protein G column, followed by dialysis against PBS.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	60868 Da
Gene Name:	BAI1 associated protein 2
Database Link:	NP_059345 Entrez Gene 108100 Mouse Entrez Gene 10458 Human Q9UQB8
Synonyms:	BAP2; FLAF3; IRSP53
Protein Families:	Druggable Genome
Protein Pathways:	Adherens junction, Regulation of actin cytoskeleton



[View online »](#)

Product images:

BAIAP2 Antibody (C-term) (Cat. #TA324869) western blot analysis in K562 cell line and mouse brain tissue lysates (35ug/lane). This demonstrates the BAIAP2 antibody detected the BAIAP2 protein (arrow).