

## Product datasheet for **TA324895**


### TPTE2 Mouse Monoclonal Antibody [Clone ID: 978CT6.1.4]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	978CT6.1.4
Applications:	WB
Recommended Dilution:	WB: 1:5000
Reactivity:	Human, Mouse
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Purified His-tagged TPIPb protein was used to produced this monoclonal antibody.
Formulation:	Mouse monoclonal antibody supplied in crude ascites with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	61112 Da
Gene Name:	transmembrane phosphoinositide 3-phosphatase and tensin homolog 2
Database Link:	<a href="#">NP_954864</a> <a href="#">Entrez Gene 93492 Human</a> <a href="#">Q6XPS3</a>
Synonyms:	TPIP
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

**Product images:**

TPIPb Antibody (C-term) (Cat. #TA324895)  
western blot analysis in 293, A431, A2058, mouse  
NIH/3T3 cell line lysates (35ug/lane). This  
demonstrates the TPIPb antibody detected the  
TPIPb protein (arrow).