


Product datasheet for **TA324906S**

USP14 Mouse Monoclonal Antibody [Clone ID: 889CT6.2.1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	889CT6.2.1
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified His-tagged USP14 protein was used to produced this monoclonal antibody.
Formulation:	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Purification:	This antibody is purified through a protein G column, followed by dialysis against PBS.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	56069 Da
Gene Name:	ubiquitin specific peptidase 14
Database Link:	NP_005142 Entrez Gene 59025 Mouse Entrez Gene 291796 Rat Entrez Gene 9097 Human P54578
Synonyms:	TGT
Protein Families:	Druggable Genome, Protease


[View online »](#)

Product images:


USP14 Antibody (N-term) (Cat. #TA324906])
western blot analysis in K562, Jurkat cell line,
mouse brain abd rat brain tissue lysates
(35ug/lane).This demonstrates the USP14
antibody detected the USP14 protein (arrow).