

## Product datasheet for **TA324943S**


### **FABP3 Mouse Monoclonal Antibody [Clone ID: 883CT19.3.3]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	883CT19.3.3
Applications:	WB
Recommended Dilution:	WB: 1:5000
Reactivity:	Mouse
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified His-tagged FABP3 protein was used to produced this monoclonal antibody.
Formulation:	Mouse monoclonal antibody supplied in crude ascites with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	14858 Da
Gene Name:	fatty acid binding protein 3
Database Link:	<a href="#">NP_004093</a> <a href="#">Entrez Gene 14077 Mouse P05413</a>
Synonyms:	FABP11; H-FABP; M-FABP; MDGI; O-FABP
Protein Pathways:	PPAR signaling pathway



[View online »](#)

**Product images:**


FAPB3 Antibody (Cat. #TA324943) western blot analysis in mouse heart tissue lysates (35ug/lane). This demonstrates the FAPB3 antibody detected the FAPB3 protein (arrow).