

Product datasheet for **TA325070S**


ASB9 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	This ASB9 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1-30 amino acids from the N-terminal region of human ASB9.
Formulation:	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Purification:	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	31858 Da
Gene Name:	ankyrin repeat and SOCS box containing 9
Database Link:	NP_001026909 Entrez Gene 140462 Human Q96DX5
Synonyms:	DKFZp564L0862; FLJ20636; MGC4954
Protein Families:	Druggable Genome



[View online »](#)

Product images:


ASB9 Antibody (N-term) (Cat. #[TA325070])
western blot analysis in A549, U251 cell line
lysates (35ug/lane). This demonstrates the ASB9
antibody detected the ASB9 protein (arrow).