


Product datasheet for **TA325152S**

ZNF434 (ZSCAN32) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	This ZNF434 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 159-187 amino acids from the Central region of human ZNF434.
Formulation:	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Purification:	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	78728 Da
Gene Name:	zinc finger and SCAN domain containing 32
Database Link:	NP_060280 Entrez Gene 54925 Human Q9NX65
Synonyms:	HCCS-5; ZNF434
Protein Families:	Transcription Factors

[View online »](#)

Product images:


ZNF434 Antibody (Center) (Cat. #TA325152]) western blot analysis in Jurkat, MCF-7 cell line lysates (35ug/lane). This demonstrates the ZNF434 antibody detected the ZNF434 protein (arrow).