

## Product datasheet for **TA325158**


### GPR58 (TAAR2) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	This TAAR2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1-30 amino acids from the N-terminal region of human TAAR2.
Formulation:	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Purification:	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	40134 Da
Gene Name:	trace amine associated receptor 2
Database Link:	<a href="#">NP_055441</a> <a href="#">Entrez Gene 9287 Human</a> <a href="#">Q9P1P5</a>
Synonyms:	GPR58; taR-2
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction



[View online »](#)

**Product images:**

TAAR2 Antibody (N-term) (Cat. #TA325158)  
western blot analysis in Jurkat cell line lysates  
(35ug/lane). This demonstrates the TAAR2  
antibody detected the TAAR2 protein (arrow).