

Product datasheet for **TA325178**


MGMT Mouse Monoclonal Antibody [Clone ID: 888CT22.2.4]

Product data:

Product Type:	Primary Antibodies
Clone Name:	888CT22.2.4
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified His-tagged MGMT protein was used to produced this monoclonal antibody.
Formulation:	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Purification:	This antibody is purified through a protein G column, followed by dialysis against PBS.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	21646 Da
Gene Name:	O-6-methylguanine-DNA methyltransferase
Database Link:	NP_002403 Entrez Gene 4255 Human P16455
Synonyms:	methylguanine-DNA methyltransferase; O-6-methylguanine-DNA methyltransferase; O6-methylguanine-DNA methyltransferase; OTTHUMP00000020741
Protein Families:	Druggable Genome



[View online »](#)

Product images:

MGMT Antibody (Cat. #TA325178) western blot analysis in MCF-7 cell line lysates (35ug/lane). This demonstrates the MGMT antibody detected the MGMT protein (arrow).