

Product datasheet for **TA325184S**

Plakophilin 4 (PKP4) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	This PKP4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1040-1068 amino acids from the C-terminal region of human PKP4.
Formulation:	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Purification:	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	131868 Da
Gene Name:	plakophilin 4
Database Link:	NP_003619 Entrez Gene 227937 Mouse Entrez Gene 8502 Human Q99569
Synonyms:	p0071



[View online »](#)

Product images:



PKP4 Antibody (C-term) (Cat. #[TA325184])
western blot analysis in 293, CEM, HL-60, K562
cell line and mouse kidney, lung tissue lysates
(35ug/lane). This demonstrates the PKP4 antibody
detected the PKP4 protein (arrow).