

Product datasheet for **TA325192S**


ZNF76 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	This ZNF76 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 388-415 amino acids from the C-terminal region of human ZNF76.
Formulation:	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Purification:	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	61831 Da
Gene Name:	zinc finger protein 76
Database Link:	NP_003418 Entrez Gene 7629 Human P36508
Synonyms:	D6S229E; Zfp523; ZNF523
Protein Families:	Transcription Factors



[View online »](#)

Product images:


ZNF76 Antibody (C-term) (Cat. #[TA325192])
western blot analysis in 293 cell line lysates
(35ug/lane). This demonstrates the ZNF76
antibody detected the ZNF76 protein (arrow).