

## **Product datasheet for TA325460**

## **Product data:**

**Product Type:** Primary Antibodies

**FGFR1 Rabbit Polyclonal Antibody** 

Applications: WE

**Recommended Dilution:** WB: 1:500-1:2000; IF/ICC: 1:100-1:500

**Reactivity:** Human, Mouse, Rat **Modifications:** Phospho-specific

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: The antiserum was produced against A synthesized peptide derived from human FGFR1

around the phosphorylation site of Tyrosine 654

Formulation: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50%

glycerol.

**Concentration:** lot specific

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific peptide.

**Conjugation:** Unconjugated

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 120 kDa

**Gene Name:** fibroblast growth factor receptor 1

Database Link: NP 001167534

Entrez Gene 14182 MouseEntrez Gene 79114 RatEntrez Gene 2260 Human

P11362

**Background:** The protein encoded by this gene is a member of the fibroblast growth factor receptor family,

where amino acid sequence is highly conserved between members and throughout

evolution. FGFR family members differ from one another in their ligand affinities and tissue

distribution.



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Synonyms: bFGF-R-1; BFGFR; CD331; CEK; ECCL; FGFBR; FGFR-1; FLG; FLT-2; FLT2; HBGFR; HH2; HRTFDS;


KAL2; N-SAM

**Protein Families:** Druggable Genome, Protein Kinase, Transmembrane

**Protein Pathways:** Adherens junction, MAPK signaling pathway, Melanoma, Pathways in cancer, Prostate cancer,

Regulation of actin cytoskeleton

## **Product images:**



Western blot analysis of FGFR1 phosphorylation expression in Insulin treated 293 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.