

Product datasheet for TA325887

STAT1 Rabbit Polyclonal Antibody

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200; IP
Reactivity:	Human, Mouse, Rat
Modifications:	Phospho-specific
Host:	Rabbit
lsotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against A synthesized peptide derived from human STAT1 around the phosphorylation site of Tyrosine 701
Formulation:	Rabbit lgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific peptide.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	84 kDa
Gene Name:	signal transducer and activator of transcription 1
Database Link:	<u>NP_009330</u> <u>Entrez Gene 20846 MouseEntrez Gene 25124 RatEntrez Gene 6772 Human</u> <u>P42224</u>
Background:	The protein encoded by this gene is a member of the STAT protein family. In response to cytokines and growth factors, STAT family members are phosphorylated by the receptor associated kinases, and then form homo- or heterodimers that translocate to the cell nucleus where they act as transcription activators.




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

STAT1 Rabbit Polyclonal Antibody – TA325887	
---	--

Synonyms:CANDF7; IMD31A; IMD31B; IMD31C; ISGF-3; STAT91Protein Families:Druggable Genome, Transcription FactorsProtein Pathways:Chemokine signaling pathway, Jak-STAT signaling pathway, Pancreatic cancer, Pathways in
cancer, Toll-like receptor signaling pathway

Product images:



Western blot analysis of STAT1 phosphorylation expression in MCF7 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US