

Product datasheet for TA325993

GRIA1 Rabbit Polyclonal Antibody

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Primary Antibodies
Applications:	IF, WB
Recommended Dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200; IF/ICC: 1:100-1:500
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
lsotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against A synthesized peptide derived from human GluR1
Formulation:	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific peptide.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	102 kDa
Gene Name:	glutamate ionotropic receptor AMPA type subunit 1
Database Link:	<u>NP_000818</u> <u>Entrez Gene 14799 MouseEntrez Gene 50592 RatEntrez Gene 2890 Human</u> <u>P42261</u>
Background:	GluR1 an integral membrane protein belonging to the glutamate-gated ion channel family. L- glutamate (Glu) acts as an excitatory neurotransmitter at many synapses in the central nervous system.
Synonyms:	GluA1; GLUH1; GLUR1; GLURA; HBGR1
Protein Families:	Druggable Genome, Ion Channels: Glutamate Receptors, Transmembrane




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

GRIA1 Rabbit Polyclonal Antibody – TA325993

Protein Pathways:

Amyotrophic lateral sclerosis (ALS), Long-term depression, Long-term potentiation, Neuroactive ligand-receptor interaction

Product images:



Western blot analysis of GluR1 expression in COS7 whole cell lysates, The lane on the left is treated with the antigen-specific peptide.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US