

## Product datasheet for **TA326059**

### Vitamin D Receptor (VDR) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000 ELISA: 1:20000
Reactivity:	Human, Mouse, Rat
Modifications:	Phospho-specific
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized phosphopeptide derived from human Vitamin D3 Receptor around the phosphorylation site of serine 51 (R-R-SP-M-K).
Formulation:	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	vitamin D (1,25- dihydroxyvitamin D3) receptor
Database Link:	<a href="#">NP_001017535</a> <a href="#">Entrez Gene 22337 Mouse</a> <a href="#">Entrez Gene 24873 Rat</a> <a href="#">Entrez Gene 7421 Human</a> <a href="#">P11473</a>
Synonyms:	NR111; PPP1R163
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors



[View online »](#)