

## Product datasheet for **TA326181**

### TrkA (NTRK1) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	IHC: 1:50~1:100 ELISA: 1:10000
Reactivity:	Human, Mouse, Rat
Modifications:	Phospho-specific
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized phosphopeptide derived from human Trk A around the phosphorylation site of tyrosine 757 (E-V-YP-A-I).
Formulation:	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	neurotrophic receptor tyrosine kinase 1
Database Link:	<a href="#">NP_001007793</a> <a href="#">Entrez Gene 18211 Mouse</a> <a href="#">Entrez Gene 59109 Rat</a> <a href="#">Entrez Gene 4914 Human</a> <a href="#">P04629</a>
Synonyms:	MTC; p140-TrkA; TRK; Trk-A; TRK1; TRKA
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane
Protein Pathways:	Apoptosis, Endocytosis, MAPK signaling pathway, Neurotrophin signaling pathway, Pathways in cancer, Thyroid cancer



[View online »](#)