

Product datasheet for **TA326234**

FRA1 (FOSL1) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000 ELISA: 1:20000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human Fra-1.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	FOS like 1, AP-1 transcription factor subunit
Database Link:	NP_005429 Entrez Gene 14283 Mouse Entrez Gene 25445 Rat Entrez Gene 8061 Human P15407
Synonyms:	FRA; fra-1; FRA1
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Wnt signaling pathway



[View online »](#)

Product images:



Western blot analysis of extracts from HT-29 cells, using Fra-1 antibody. The lane on the right is treated with the synthesized peptide.