

Product datasheet for **TA326375**

Hsp60 (HSPD1) Rabbit Polyclonal Antibody [Clone ID: N/A]

Product data:

Product Type:	Primary Antibodies
Clone Name:	N/A
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human, Mouse, Rat, Bovine, Canine, Chicken, Rabbit
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Human Hsp60 produced through recombinant DNA methods in E.coli
Formulation:	PBS with 50% glycerol and 0.09% sodium azide
Concentration:	lot specific
Purification:	Protein A Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	heat shock protein family D (Hsp60) member 1
Database Link:	NP_002147 Entrez Gene 15510 Mouse Entrez Gene 63868 Rat Entrez Gene 3329 Human P10809



[View online »](#)

Background:

In both prokaryotic and eukaryotic cells, the misfolding and aggregation of proteins during biogenesis and under conditions of cellular stress are prevented by molecular chaperones. Members of the HSP60 family of heat shock proteins are some of the best characterized chaperones. Hsp60, also known as Cpn60 or GroEL, is an abundant protein synthesized constitutively in the cell that is induced to a higher concentration after brief cell shock. It is present in many species and exhibits a remarkable sequence homology among various counterparts in bacteria, plants, and mammals with more than half of the residues identical between bacterial and mammalian Hsp60. Whereas mammalian Hsp60 is localized within the mitochondria, plant Hsp60, or otherwise known as Rubisco-binding protein, is located in plant chloroplasts. It has been indicated that these proteins carry out a very important biological function due to the fact that Hsp60 is present in so many different species. The common characteristics of the Hsp60s from the divergent species are i) high abundance, ii) induction with environmental stress such as heat shock, iii) homo-oligomeric structures of either 7 or 14 subunits which reversibly dissociate in the presence of Mg^{2+} and ATP, iv) ATPase activity and v) a role in folding and assembly of oligomeric protein structures. These similarities are supported by recent studies where the single-ring human mitochondrial homolog, Hsp60 with its co-chaperonin, Hsp10 were expressed in a *E. coli* strain, engineered so that the *groE* operon is under strict regulatory control. This study has demonstrated that expression of Hsp60-Hsp10 was able to carry out all essential *in vivo* functions of GroEL and its co-chaperonin, GroES. Hsp60 has however been linked to a number of autoimmune diseases, as well as Alzheimers, coronary artery diseases, MS, and diabetes.

Synonyms:

CPN60; GROEL; HLD4; HSP-60; HSP60; HSP65; HuCHA60; SPG13

Note:


Recognizes a ~60kDa protein, corresponding to the molecular mass of hsp60 on SDS-PAGE blots.

Protein Families:

Druggable Genome, Stem cell - Pluripotency

Protein Pathways:

RNA degradation, Type I diabetes mellitus

Product images:

Western blot analysis of Hsp60 using a 1:1000 dilution of the antibody