

Product datasheet for TA327542

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

H3FT (HIST3H3) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: ChIP, ChIP-seq, ICC/IF, IP, WB

Recommended Dilution: WB 1:500 - 1:2000;IF 1:50 - 1:200;IP 1:50 - 1:200;ChIP 1:20 - 1:100

Reactivity: Human, Mouse, Rat, Other (Wide Range)

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: A synthetic methylated peptide corresponding to residues surrounding R2 of human histone

Н3

Formulation: Store at -20C or -80C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50%

glycerol, pH7.3

Concentration: lot specific

Purification: Affinity purification

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 16 kDa

Gene Name: histone cluster 3, H3

Database Link: NP 003484

Entrez Gene 691496 RatEntrez Gene 8290 Human

Q16695






Background:

Modulation of chromatin structure plays an important role in the regulation of transcription in eukaryotes. The nucleosome, made up of DNA wound around eight core histone proteins (two each of H2A, H2B, H3, and H4), is the primary building block of chromatin. The aminoterminal tails of core histones undergo various post-translational modifications, including acetylation, phosphorylation, methylation, and ubiquitination. These modifications occur in response to various stimuli and have a direct effect on the accessibility of chromatin to transcription factors and, therefore, gene expression. In most species, histone H2B is primarily acetylated at Lys5, 12, 15, and 20. Histone H3 is primarily acetylated at Lys9, 14, 18, 23, 27, and 56. Acetylation of H3 at Lys9 appears to have a dominant role in histone deposition and chromatin assembly in some organisms. Phosphorylation at Ser10, Ser28, and Thr11 of histone H3 is tightly correlated with chromosome condensation during both mitosis and meiosis. Phosphorylation at Thr3 of histone H3 is highly conserved among many species and is catalyzed by the kinase haspin. Immunostaining with phospho-specific antibodies in mammalian cells reveals mitotic phosphorylation at Thr3 of H3 in prophase and its dephosphorylation during anaphase.

Synonyms: g; H3; H3.4; H3FT; H3t


Protein Pathways: Systemic lupus erythematosus

Product images:



Western blot analysis of extracts of various cell lines, using H3R2me2a antibody.





Immunofluorescence analysis of 293T cell using H3R2me2a antibody. Blue: DAPI for nuclear staining.