

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA327652**

# CD71 (TFRC) Mouse Monoclonal Antibody [Clone ID: MRQ-48]

# **Product data:**

Product Type:	Primary Antibodies
Clone Name:	MRQ-48
Applications:	IHC
Recommended Dilution:	IHC: 1:100 - 1:500
Reactivity:	Human
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Formulation:	This antibody is supplied as cell culture supernatant diluted in tris buffered saline, pH 7.3-7.7, with 1% BSA and <0.1% sodium azide.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	transferrin receptor
Database Link:	<u>NP 003225</u> <u>Entrez Gene 7037 Human</u> <u>P02786</u>
Synonyms:	CD71; IMD46; p90; T9; TFR; TFR1; TR; TRFR



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

### **GRIGENE** CD71 (TFRC) Mouse Monoclonal Antibody [Clone ID: MRQ-48] – TA327652

The transferrin receptor (CD71) is most highly expressed on placental syncytiotrophoblasts, Note: myocytes, basal keratinocytes, hepatocytes, endocrine pancreas, spermatocytes, and erythroid precursors. The level of transferrin receptor expression is highest in early erythroid precursors through the intermediate normoblast phase, after which expression decreases through the reticulocyte phase. The maturation of erythrocytes results in loss of transferrin receptor expression, in concert with down-regulation of the machinery for hemoglobin synthesis. The high level of transferrin receptor within erythroid precursors makes anti-CD71 an excellent marker for evaluation of erythroid precursors within bone marrow biopsy specimens and shows the following features: 1) distinct membranous and cytoplasmic staining pattern, which is easily recognized in bone marrow biopsy; 2) restriction to erythroid lineage within bone marrow biopsy specimens; 3) CD71 expression decreases with the maturation of erythrocytes, with the highest level seen in early forms and the lowest level in late normoblast stage, and most importantly; 4) mature erythrocytes do not express CD71, which facilitates bone marrow analyses. Anti-CD71 is useful in identifying erythroid precursors with very little interference from mature erythrocytes and also in the determination of erythroid leukemia, benign erythroid proliferative disorders, and myelodysplastic syndrome, although further studies are needed for making a definitive diagnosis of myelodysplastic syndrome.

# Protein Families:Druggable Genome, ES Cell Differentiation/IPS, Protease, Secreted Protein, TransmembraneProtein Pathways:Endocytosis, Hematopoietic cell lineage

### **Product images:**



Immunohistochemistry staining of Paraffin Bone marrow tissue by CD71 antibody (dilution: 1:100 -1:500; visualization of staining: Cytoplasmic, membranous)

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US