

## Product datasheet for **TA327763**


### Fibrinogen beta chain (FGB) Goat Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 0.01-0.03 ug/mL
Reactivity:	Human (Expected from sequence similarity: Pig)
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	internal region (EKHQLYIDETVNSN)
Formulation:	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Concentration:	lot specific
Purification:	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	fibrinogen beta chain
Database Link:	<a href="#">NP_005132</a> <a href="#">Entrez Gene 2244 Human</a> <a href="#">P02675</a>
Synonyms:	HEL-S-78p
Note:	Approx 60kDa band observed in Human Liver lysates (calculated MW of 55.9kDa according to NP_005132.2). Recommended concentration: 0.01-0.03µg/ml.
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Complement and coagulation cascades



[View online »](#)

**Product images:**

Anti-FGB antibody (0.01ug/ml) staining of Human Liver lysate (35ug protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.