

Product datasheet for **TA328112S**

ATP5PB Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	This ATP5F1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 161-195 amino acids from the Central region of human ATP5F1.
Formulation:	PBS with 0.09% (W/V) sodium azide
Concentration:	lot specific
Purification:	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ATP synthase, H ⁺ transporting, mitochondrial Fo complex subunit B1
Database Link:	NP_001679 Entrez Gene 11950 Mouse Entrez Gene 171375 Rat Entrez Gene 515 Human P24539
Background:	Mitochondrial membrane ATP synthase (F(1)F(0) ATP synthase or Complex V) produces ATP from ADP in the presence of a proton gradient across the membrane which is generated by electron transport complexes of the respiratory chain. F-type ATPases consist of two structural domains, F(1) -containing the extramembraneous catalytic core, and F(0) - containing the membrane proton channel, linked together by a central stalk and a peripheral stalk. During catalysis, ATP synthesis in the catalytic domain of F(1) is coupled via a rotary mechanism of the central stalk subunits to proton translocation. Part of the complex F(0) domain and the peripheric stalk, which acts as a stator to hold the catalytic alpha(3)beta(3) subcomplex and subunit a/ATP6 static relative to the rotary elements.


[View online »](#)

Synonyms: PIG47

Protein Pathways: Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease

Product images:



ATP5F1 Antibody (Center) (Cat. #TA328112J)
western blot analysis in HepG2, HL-60 cell line,
mouse kidney and rat heart tissue lysates
(35ug/lane). This demonstrates the ATP5F1
antibody detected the ATP5F1 protein (arrow).