

Product datasheet for **TA328130**


PPHLN1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	This PPHLN1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 45-66 amino acids from the N-terminal region of human PPHLN1.
Formulation:	PBS with 0.09% (W/V) sodium azide
Concentration:	lot specific
Purification:	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	periphilin 1
Database Link:	NP_958923 Entrez Gene 51535 Human Q8NEY8
Background:	Involved in epithelial differentiation and contributes to epidermal integrity and barrier formation.
Synonyms:	CR; HSPC206; HSPC232



[View online »](#)

Product images:

PPHLN1 Antibody (N-term) (Cat. #TA328130) western blot analysis in A549 cell line lysates (35ug/lane). This demonstrates the PPHLN1 antibody detected the PPHLN1 protein (arrow).