


Product datasheet for **TA328189S**

GGN Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:1000
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	This GGN antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 139-163 amino acids from the Central region of human GGN.
Formulation:	PBS with 0.09% (W/V) sodium azide
Concentration:	lot specific
Purification:	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	gametogenetin
Database Link:	NP_872283 Entrez Gene 199720 Human Q86UU5
Background:	May be involved in spermatogenesis (By similarity).
Synonyms:	FLJ35713; gametogenetin; MGC33369

[View online »](#)

Product images:


GGN Antibody (Center) (Cat. #[TA328189])
western blot analysis in A549, LNCaP cell line
lysates (35ug/lane). This demonstrates the GGN
antibody detected the GGN protein (arrow).