

Product datasheet for **TA328219**

IL7 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	WB: 0.1-0.2ug/mL, ELISA: 0.25-2ug/mL
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	E.coli derived Recombinant Human IL-7
Formulation:	A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.
Purification:	Produced from sera of rabbits pre-immunized with highly pure recombinant Human IL-7. Anti-Human IL-7 specific antibody was purified by affinity chromatography employing immobilized Human IL-7 matrix.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	interleukin 7
Database Link:	NP_000871 Entrez Gene 3574 Human P13232
Synonyms:	IL-7
Note:	Neutralization: To yield one-half maximal inhibition [ND50] of the biological activity of Human IL-7 (1.0ng/ml), a concentration of 0.04 - 0.06 <=ug/ml of this antibody is required.
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling pathway




[View online »](#)

Product images:



To detect Human IL-7 by sandwich ELISA (using 100 ul/well antibody solution) a concentration of 0.5 - 2.0 ug/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Human IL-7 ([TA328218]) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant Human IL-7.



To detect Human IL-7 by sandwich ELISA (using 100 ul/well antibody solution) a concentration of 0.5 - 2.0 ug/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Human IL-7 ([TA328218]) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant Human IL-7.