

## Product datasheet for **TA328280**


### CNTF Rabbit Polyclonal Antibody

#### Product data:


Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	WB: 0.1-0.2ug/mL, ELISA: 0.25-2ug/mL, IHC: 0.25ug/mL - 1mg/mL
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	E.coli derived Recombinant Human CNTF
Formulation:	A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.
Purification:	Produced from sera of rabbits pre-immunized with highly pure (>98%) recombinant hCNTF. Anti-Human CNTF specific antibody was purified by affinity chromatography employing immobilized hCNTF matrix.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ciliary neurotrophic factor
Database Link:	<a href="#">NP_000605</a> <a href="#">Entrez Gene 1270 Human P26441</a>
Synonyms:	HCNTF
Note:	Neutralization: To yield one-half maximal inhibition [ND50] of the biological activity of Human CNTF (1.0 ug/ml) a concentration of 7.0 - 11.0 ug/ml of this antibody is required.
Protein Families:	Druggable Genome
Protein Pathways:	Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway




[View online »](#)

**Product images:**


To detect hCNTF by sandwich ELISA (using 100 ul/well antibody solution) a concentration of 0.5 - 2.0 ug/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Human CNTF ([TA328279]) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hCNTF.



To detect hCNTF by sandwich ELISA (using 100 ul/well antibody solution) a concentration of 0.5 - 2.0 ug/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Human CNTF ([TA328279]) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hCNTF.



To detect hCNTF by sandwich ELISA (using 100 ul/well antibody solution) a concentration of 0.5 - 2.0 ug/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Human CNTF ([TA328279]) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hCNTF.