

Product datasheet for **TA328296**


Prokineticin 1 (PROK1) Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	WB: 0.1-0.2ug/mL, ELISA: 0.25-2ug/mL
Reactivity:	Human
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	E.coli derived Recombinant Human EG-VEGF
Formulation:	A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.
Purification:	Produced from sera of goats pre-immunized with highly pure (>98%) recombinant hEG-VEGF. Anti-Human EG-VEGF specific antibody was purified by affinity chromatography employing immobilized hVEGF matrix.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	prokineticin 1
Database Link:	NP_115790 Entrez Gene 84432 Human P58294
Synonyms:	EGVEGF; PK1; PRK1
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

Product images:

To detect hEG-VEGF by sandwich ELISA (using 100 ul/well antibody solution) a concentration of 0.5 - 2.0 ug/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Human EG-VEGF ([TA328295]) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hEG-VEGF.