

Product datasheet for **TA328473**


Visfatin (NAMPT) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	WB: 0.1-0.2ug/mL, ELISA: 0.25-2ug/mL, IHC: 7.5-15ng/mL
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	E.coli derived Recombinant Human Visfatin
Formulation:	A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.
Purification:	Produced from sera of rabbits immunized with highly pure recombinant Human Visfatin. Anti-Human Visfatin specific antibody was purified by affinity chromatography employing an immobilized Human Visfatin matrix.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	nicotinamide phosphoribosyltransferase
Database Link:	NP_005737 Entrez Gene 10135 Human P43490
Synonyms:	1110035O14Rik; PBEF; PBEF1; VF; VISFATIN
Protein Families:	Druggable Genome
Protein Pathways:	Nicotinate and nicotinamide metabolism




[View online »](#)

Product images:


To detect Human Visfatin by sandwich ELISA (using 100ul/well antibody solution) a concentration of 0.5 - 2.0 ug/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Human Visfatin ([TA328472]) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant Human Visfatin.



To detect Human Visfatin by sandwich ELISA (using 100ul/well antibody solution) a concentration of 0.5 - 2.0 ug/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Human Visfatin ([TA328472]) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant Human Visfatin.



To detect Human Visfatin by sandwich ELISA (using 100ul/well antibody solution) a concentration of 0.5 - 2.0 ug/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Human Visfatin ([TA328472]) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant Human Visfatin.