

## Product datasheet for **TA328483**

### IL3 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	WB: 0.1-0.2ug/mL, ELISA: 0.25-2ug/mL
Reactivity:	Murine
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	E.coli derived Recombinant Murine IL-3
Formulation:	A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.
Purification:	Produced from sera of rabbits pre-immunized with highly pure (>98%) recombinant mIL-3. Anti-Murine IL-3 specific antibody was purified by affinity chromatography employing immobilized mIL-3 matrix.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	interleukin 3
Database Link:	<a href="#">NP_034686</a> <a href="#">Entrez Gene 16187 Mouse</a> <a href="#">P01586</a>
Synonyms:	IL-3; MCGF; MGC79398; MGC79399; MULTI-CSF
Note:	Neutralization: To yield one-half maximal inhibition [ND50] of the biological activity of mIL-3 (0.5 ng/ml), a concentration of 0.15 - 0.20 ≤ug/ml of this antibody is required.




[View online »](#)

Product images:



To detect mIL-3 by sandwich ELISA (using 100 ul/well antibody solution) a concentration of 0.5 - 2.0 ug/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Murine IL-3 ([TA328482]) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant mIL-3.



To detect mIL-3 by sandwich ELISA (using 100 ul/well antibody solution) a concentration of 0.5 - 2.0 ug/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Murine IL-3 ([TA328482]) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant mIL-3.