

Product datasheet for **TA328487**

IL5 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	WB: 0.1-0.2ug/mL, ELISA: 0.25-2ug/mL
Reactivity:	Murine
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	E.coli derived Recombinant Murine IL-5
Formulation:	A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.
Purification:	Produced from sera of rabbits immunized with highly pure recombinant Murine IL-5. Anti-Murine IL-5 specific antibody was purified by affinity chromatography employing an immobilized Murine IL-5 matrix.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	interleukin 5
Database Link:	NP_034688 Entrez Gene 16191 Mouse P04401
Synonyms:	EDF; IL-5; interleukin-5; TRF




[View online »](#)

Product images:



To detect Murine IL-5 by sandwich ELISA (using 100ul/well antibody solution) a concentration of 0.5 - 2.0 ug/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Murine IL-5 ([TA328486]) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant Murine IL-5.



To detect Murine IL-5 by sandwich ELISA (using 100ul/well antibody solution) a concentration of 0.5 - 2.0 ug/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Murine IL-5 ([TA328486]) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant Murine IL-5.