

Product datasheet for **TA328511**

Ccl11 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	WB: 0.1-0.2ug/mL, ELISA: 0.25-2ug/mL
Reactivity:	Murine
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	E.coli derived Recombinant Murine Eotaxin (CCL11)
Formulation:	A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.
Purification:	Produced from sera of rabbits pre-immunized with highly pure (>98%) recombinant mEotaxin. Anti-Murine Eotaxin specific antibody was purified by affinity chromatography employing immobilized mEotaxin matrix.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	chemokine (C-C motif) ligand 11
Database Link:	NP_035460 Entrez Gene 20292 Mouse P48298
Synonyms:	eotaxin; eotaxin-1; MGC22554; SCYA11
Note:	Neutralization: To yield one-half maximal inhibition [ND50] of the biological activity of mEotaxin (100.00 ng/ml), a concentration of 5.0-8.0 ≤ug/ml of this antibody is required.




[View online »](#)

Product images:



To detect mEotaxin by sandwich ELISA (using 100 ul/well antibody solution) a concentration of 0.5 - 2.0 ug/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Murine Eotaxin ([TA328510]) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant mEotaxin.



To detect mEotaxin by sandwich ELISA (using 100 ul/well antibody solution) a concentration of 0.5 - 2.0 ug/ml of this antibody is required. This antigen affinity purified antibody, in conjunction with Biotinylated Anti-Murine Eotaxin ([TA328510]) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant mEotaxin.