

Product datasheet for **TA328592**

Ccl5 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	WB: 0.1-0.2ug/mL, ELISA: 0.25-2ug/mL
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	E.coli derived Recombinant Rat RANTES (CCL5)
Formulation:	A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.
Purification:	Produced from sera of rabbits pre-immunized with highly pure (>98%) recombinant Rat RANTES. Anti-Rat RANTES specific antibody was purified by affinity chromatography and then biotinylated.
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	C-C motif chemokine ligand 5
Database Link:	NP_112378 Entrez Gene 81780 Rat
Synonyms:	D17S136E; EoCP; MGC17164; OTTHUMP00000197106; RANTES; SCYA5; SIS-delta; SISd; TCP228




[View online »](#)

Product images:



To detect Rat RANTES by direct ELISA (using 100 ul/well antibody solution) a concentration of 0.25-1.0 ug/ml of this antibody is required. This biotinylated polyclonal antibody, in conjunction with compatible secondary reagents, allows the detection of at least 0.2-0.4 ng/well of recombinant Rat RANTES.



To detect Rat RANTES by sandwich ELISA (using 100 ul/well antibody solution) a concentration of 0.25-1.0 ug/ml of this antibody is required. This biotinylated polyclonal antibody, in conjunction with Polyclonal Anti-Rat RANTES ([TA328593]) as a capture antibody, allows the detection of at least 0.2-0.4 ng/well of recombinant Rat RANTES.