

Product datasheet for **TA329140**

p53 (TP53) Rabbit Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB, IHC, ChIP
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-TP53 antibody: synthetic peptide directed towards the N terminal of human TP53. Synthetic peptide located within the following region: EEPQSDPSVEPPLSQETFSDLWKLLPENNVLSPLPSQAMDDLMLSPDDIE
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	44 kDa
Gene Name:	tumor protein p53
Database Link:	NP_000537 Entrez Gene 7157 Human P04637




[View online »](#)

Background:	TP53 acts as a tumor suppressor in many tumor types; induces growth arrest or apoptosis depending on the physiological circumstances and cell type. Involved in cell cycle regulation as a trans-activator that acts to negatively regulate cell division by controlling a set of genes required for this process. This gene encodes tumor protein p53, which responds to diverse cellular stresses to regulate target genes that induce cell cycle arrest, apoptosis, senescence, DNA repair, or changes in metabolism. p53 protein is expressed at low level in normal cells and at a high level in a variety of transformed cell lines, where it's believed to contribute to transformation and malignancy. p53 is a DNA-binding protein containing transcription activation, DNA-binding, and oligomerization domains. It is postulated to bind to a p53-binding site and activate expression of downstream genes that inhibit growth and/or invasion, and thus function as a tumor suppressor. Mutants of p53 that frequently occur in a number of different human cancers fail to bind the consensus DNA binding site, and hence cause the loss of tumor suppressor activity. Alterations of this gene occur not only as somatic mutations in human malignancies, but also as germline mutations in some cancer-prone families with Li-Fraumeni syndrome. Multiple p53 variants due to alternative promoters and multiple alternative splicing have been found. These variants encode distinct isoforms, which can regulate p53 transcriptional activity.
Synonyms:	BCC7; LFS1; P53; TRP53
Note:	Immunogen sequence homology: Pig: 100%; Human: 100%; Sheep: 100%; Bovine: 100%; Rabbit: 100%; Rat: 90%; Horse: 90%; Guinea pig: 90%; Dog: 80%
Protein Families:	Druggable Genome, Stem cell - Pluripotency, Transcription Factors
Protein Pathways:	Amyotrophic lateral sclerosis (ALS), Apoptosis, Basal cell carcinoma, Bladder cancer, Cell cycle, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, Glioma, Huntington's disease, MAPK signaling pathway, Melanoma, Neurotrophin signaling pathway, Non-small cell lung cancer, p53 signaling pathway, Pancreatic cancer, Pathways in cancer, Prostate cancer, Small cell lung cancer, Thyroid cancer, Wnt signaling pathway


Product images:




WB Suggested Anti-TP53 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1:1562500; Positive Control: 293T cell lysate. TP53 is strongly supported by BioGPS gene expression data to be expressed in HEK293T




Host: Rabbit; Target Name: TP53; Sample Tissue: Human Fetal Liver; Antibody Dilution: 1.0ug/ml




Host: Rabbit; Target Name: TP53; Sample Tissue: Human Adult Placenta; Antibody Dilution: 1.0ug/ml




Host: Rabbit; Target Name: TP53; Sample Tissue: Human 721_B; Antibody Dilution: 1.0ug/ml; TP53 is strongly supported by BioGPS gene expression data to be expressed in Human 721_B cells



Immunohistochemistry with U20S/ PLKO cells tissue



Immunohistochemistry with U2OS/ PLKO cells tissue



Sample Type : U2OS/ PLKO cells

Figure 1. Binding of p53-specific antibodies to the p21 promoter.

Application: ChIP Assay
 Application data in forum

Submitted by:
 Nick Barlev
 Department of Biochemistry
 University of Leicester