

Product datasheet for **TA329227**

ZNF187 (ZSCAN26) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF187 antibody: synthetic peptide directed towards the n terminal of human ZNF187. Synthetic peptide located within the following region: IILPKELQARVQEHHPEsREDVVVLEDLQLDLGETGQQDPDQPKKQKIL
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	55 kDa
Gene Name:	zinc finger and SCAN domain containing 26
Database Link:	NP_001104509 Entrez Gene 7741 Human Q16670
Background:	ZNF187 may be involved in transcriptional regulation.
Synonyms:	SRE-ZBP; ZNF187
Note:	Immunogen sequence homology: Human: 100%; Pig: 85%
Protein Families:	Transcription Factors



[View online »](#)

Product images:



WB Suggested Anti-ZNF187 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:62500; Positive Control:
SH-SY5Y cell lysate