

Product datasheet for **TA329370**


ZNF449 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-ZNF449 antibody: synthetic peptide directed towards the middle region of human ZNF449. Synthetic peptide located within the following region: LENREEPWVKELQDSKEMKQLLDSKIGFEIGIENEEDTSKQKKMETMYPF
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	60 kDa
Gene Name:	zinc finger protein 449
Database Link:	NP_689908 Entrez Gene 203523 Human Q6P9G9
Background:	ZNF449 belongs to the krueppel C2H2-type zinc-finger protein family. It contains 7 C2H2-type zinc fingers and 1 SCAN box domain. ZNF449 may be involved in transcriptional regulation.
Synonyms:	ZSCAN19
Note:	Immunogen sequence homology: Dog: 100%; Pig: 100%; Human: 100%; Bovine: 100%; Horse: 93%; Mouse: 93%; Rabbit: 93%; Rat: 86%; Guinea pig: 77%
Protein Families:	Transcription Factors



[View online »](#)

Product images:

WB Suggested Anti-ZNF449 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:312500; Positive
Control: Transfected 293T