

Product datasheet for **TA329404**

Foxs1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-Fkhl18 antibody: synthetic peptide corresponding to a region of Mouse. Synthetic peptide located within the following region: MQKQPSPESLAPSAEPTKPPYSYIALIAMAIQSSPGQRATLSGIYRYIMG
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	35 kDa
Gene Name:	forkhead box S1
Database Link:	NP_034356 Entrez Gene 14239 Mouse Q61574
Background:	Fkhl18 is a transcriptional repressor that suppresses transcription from the FASLG, FOXO3 and FOXO4 promoters. Fkhl18 may have a role in the organization of the testicular vasculature.
Synonyms:	FKHL18; FREAC-10; FREAC10; MGC4544
Note:	Immunogen sequence homology: Pig: 100%; Human: 100%; Guinea pig: 100%; Rat: 92%; Horse: 85%; Mouse: 85%; Bovine: 85%; Rabbit: 85%



[View online »](#)

Product images:



WB Suggested Anti-Fkhl18 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1:312500; Positive Control: Mouse Liver