

Product datasheet for **TA329489**


Pla2g12a Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-Pla2g12a antibody is: synthetic peptide directed towards the middle region of Rat Pla2g12a. Synthetic peptide located within the following region: IRNGIHKIDTYLNAALDLLGGEDGLCQYKCSGSKPAPRYGYKPSPPNGC
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	12 kDa
Gene Name:	phospholipase A2, group XIA
Database Link:	NP_001102035 Entrez Gene 362039 Rat
Background:	The function of this protein remains unknown.
Synonyms:	GXII; PLA2G12; ROSSY; sPLA2-XII
Note:	Immunogen sequence homology: Dog: 100%; Pig: 100%; Rat: 100%; Human: 100%; Mouse: 100%; Rabbit: 100%; Guinea pig: 100%; Bovine: 93%



[View online »](#)

Product images:

Host: Rabbit; Target Name: Pla2g12a; Sample
Tissue: Rat Spleen lysates; Antibody Dilution:
1.0ug/ml