

## Product datasheet for **TA329784**

### **C10orf30 (BEND7) Rabbit Polyclonal Antibody**


#### **Product data:**

|                         |   |
|-------------------------|---|
| Product Type:           | Primary Antibodies  |
| Applications:           | WB  |
| Recommended Dilution:   | WB  |
| Reactivity:             | Human   |
| Host:                   | Rabbit  |
| Isotype:                | IgG   |
| Clonality:              | Polyclonal  |
| Immunogen:              | The immunogen for anti-C10orf30 antibody: synthetic peptide directed towards the N terminal of human C10orf30. Synthetic peptide located within the following region: AGSNCTCNCQSTLQAILQELKTMRKLMQIQAVGTQNRQQPPISLICSQR |
| Formulation:            | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.<br><i>Note that this product is shipped as lyophilized powder to China customers.</i>                                 |
| Purification:           | Affinity Purified   |
| Conjugation:            | Unconjugated  |
| Storage:                | Store at -20°C as received.   |
| Stability:              | Stable for 12 months from date of receipt.  |
| Predicted Protein Size: | 51 kDa  |
| Gene Name:              | BEN domain containing 7   |
| Database Link:          | <a href="#">NP_689964</a><br><a href="#">Entrez Gene 222389 Human</a><br><a href="#">Q8N7W2</a>   |
| Background:             | The function of this protein remains unknown.   |
| Synonyms:               | C10orf30  |
| Note:                   | Immunogen sequence homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%   |



[View online »](#)

## Product images:



WB Suggested Anti-C10orf30 Antibody Titration:  
0.2-1 ug/ml; Positive Control: HT1080 cell lysate