

Product datasheet for TA329946

Dmrtc2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-DMRTC2 antibody: synthetic peptide corresponding to a region of Mouse. Synthetic peptide located within the following region: MDPSETAALHHCSADSSPADEARVPQSTELIPRRPVRSRPTCARCRNHGV
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	41 kDa
Gene Name:	doublesex and mab-3 related transcription factor like family C2
Database Link:	NP_082008 Entrez Gene 71241 Mouse Q8CGW9
Background:	Dmrtc2 may be involved in sexual development
Synonyms:	DMRTC2
Note:	Immunogen sequence homology: Rat: 100%; Mouse: 100%



[View online »](#)

Product images:



WB Suggested Anti-DMRTC2 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1:62500; Positive Control:
Human Testis