

Product datasheet for **TA329957**

Pbx4 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-Pbx4 antibody: synthetic peptide corresponding to a region of Mouse. Synthetic peptide located within the following region: MAAPLRVPPQPAPRRRLPTTAPLGHDTSDVLQQIMAITDQSLDEAQARKH
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	42 kDa
Gene Name:	pre B cell leukemia homeobox 4
Database Link:	NP_001020125 Entrez Gene 80720 Mouse Q99NE9
Background:	The function of this protein remains unknown.
Synonyms:	PBX4
Note:	Immunogen sequence homology: Mouse: 100%; Rat: 92%



[View online »](#)

Product images:



WB Suggested Anti-Pbx4 Antibody Titration: 0.2-1 ug/ml; Positive Control: Mouse Small Intestine