

## Product datasheet for TA330463

### Unkl Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-Unkl antibody is: synthetic peptide directed towards the middle region of Mouse Unkl. Synthetic peptide located within the following region: IRQWEESWQQVKQACDAWQREAQEAKERARVADSDRQLALQRKEEVEAKV
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	26 kDa
Gene Name:	unkempt family like zinc finger
Database Link:	<a href="#">NP_083065</a> <a href="#">Entrez Gene 74154 Mouse</a> <a href="#">Q5FWH2</a>
Background:	The function of this protein remains unknown.
Synonyms:	C16orf28; FLJ12623; FLJ23360; KIAA0734; MGC5179; OTTHUMP00000197302; unkempt-like; ZC3H5L; ZC3HDC5L
Note:	Immunogen sequence homology: Pig: 100%; Rat: 100%; Human: 100%; Mouse: 100%; Guinea pig: 100%



[View online »](#)

Product images:



Host: Rabbit; Target Name: Unkl; Sample Tissue: Mouse Stomach lysates; Antibody Dilution: 1.0ug/ml