

Product datasheet for **TA330674**


TMEM143 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-TMEM143 antibody is: synthetic peptide directed towards the C-terminal region of Human TMEM143. Synthetic peptide located within the following region: VTFNGTRALAHQLQALTPSMGLYPPPGFKLDPVAPITSEPPQATPSSNIS
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	50 kDa
Gene Name:	transmembrane protein 143
Database Link:	NP_060743 Entrez Gene 55260 Human Q96AN5
Background:	The function of this protein remains unknown.
Synonyms:	FLJ10922
Note:	Human: 100%
Protein Families:	Transmembrane


[View online »](#)

Product images:



Host: Rabbit; Target Name: TMEM143; Sample
Tissue: Jurkat Whole Cell lysates; Antibody
Dilution: 1.0ug/ml