

Product datasheet for **TA330706**

Pdzd4 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-Pdzd4 antibody is: synthetic peptide directed towards the C-terminal region of Mouse Pdzd4. Synthetic peptide located within the following region: SEMKMGRYWSKEERKQHLIRAREQRKRREFMMQSRLECLREQQNGDSKPE
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	86 kDa
Gene Name:	PDZ domain containing 4
Database Link:	NP_001025039 Entrez Gene 245469 Mouse Q9QY39
Background:	The function of this protein remains unknown.
Synonyms:	FLJ34125; KIAA1444; LU1; PDZK4; PDZRN4L
Note:	Rat: 100%; Human: 100%; Mouse: 100%; Dog: 93%; Pig: 93%; Horse: 93%; Bovine: 93%; Guinea pig: 93%



[View online »](#)

Product images:



Host: Rabbit; Target Name: Pdzd4; Sample
Tissue: Mouse Heart lysates; Antibody Dilution:
1.0ug/ml