

## Product datasheet for **TA330885**

### PRAMEF19 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-PRAMEF19 antibody is: synthetic peptide directed towards the N-terminal region of Human PRAMEF19. Synthetic peptide located within the following region: EKQPLKVFMDVCLKEKSVDELDLFFSGWVQHRRRSVHLCCTKVVNYSMNI
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	45 kDa
Gene Name:	PRAME family member 19
Database Link:	<a href="#">NP_001093260</a> <a href="#">Entrez Gene 645414 Human</a> <a href="#">Q5SWL8</a>
Background:	The function of this protein remains unknown.
Synonyms:	PRAMEF18
Note:	Human: 100%



[View online »](#)

## Product images:



Host: Rabbit  
Target Name: PRAMEF19  
Sample Tissue: 721\_B Cell Lysate  
Antibody Dilution: 1.0µg/ml

Host: Rabbit; Target Name: PRAMEF19; Sample  
Tissue: 721\_B Whole Cell lysates; Antibody  
Dilution: 1.0ug/ml