

Product datasheet for **TA330930**


PXT1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for anti-PXT1 antibody: synthetic peptide directed towards the middle region of human PXT1. Synthetic peptide located within the following region: MQLRHIGDNIDHRMVREDLQQDGRDALDHFVFFFRRVQVLLHFFWNNHL
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Purification:	Affinity Purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	6 kDa
Gene Name:	peroxisomal, testis specific 1
Database Link:	NP_694535 Entrez Gene 222659 Human Q8NFPO
Background:	The exact function of PXT1 remains unknown.
Synonyms:	STEPP
Note:	Human: 100%; Horse: 92%



[View online »](#)

Product images:

WB Suggested Anti-PXT1 Antibody Titration: 0.2-1 ug/ml; ELISA Titer: 1:312500; Positive Control: Jurkat cell lysate