

Product datasheet for **TA331504**

ZRANB2 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-ZRANB2 Antibody: synthetic peptide directed towards the N terminal of human ZRANB2. Synthetic peptide located within the following region: KTCSNVNWARRSECNMCNTPKYAKLEERTGYGGGFNERENVEYIEREESD
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	37 kDa
Gene Name:	zinc finger RANBP2-type containing 2
Database Link:	NP_976225 Entrez Gene 9406 Human O95218
Background:	ZRANB2 is a splice factor required for alternative splicing of SFRS10/TRA2B transcripts. It may interfere with constitutive 5'-splice site selection.
Synonyms:	ZIS; ZIS1; ZIS2; ZNF265
Note:	Immunogen sequence homology: Dog: 100%; Pig: 100%; Rat: 100%; Horse: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%; Yeast: 90%; Zebrafish: 79%
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

Product images:



WB Suggested Anti-ZRANB2 Antibody Titration:
0.2-1 ug/ml; ELISA Titer: 1: 312500; Positive
Control: MCF7 cell lysate