

Product datasheet for **TA331621**

ZNF449 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen for Anti-ZNF449 Antibody: synthetic peptide directed towards the N terminal of human ZNF449. Synthetic peptide located within the following region: VLEQFLTILPTEIETWVREHCPENRERVSLIEDLQRELEIPEQQVDMHD
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	60 kDa
Gene Name:	zinc finger protein 449
Database Link:	NP_689908 Entrez Gene 203523 Human Q6P9G9
Background:	ZNF449 belongs to the krueppel C2H2-type zinc-finger protein family. It contains 7 C2H2-type zinc fingers and 1 SCAN box domain. ZNF449 may be involved in transcriptional regulation.
Synonyms:	ZSCAN19
Note:	Immunogen sequence homology: Dog: 100%; Pig: 100%; Rat: 100%; Human: 100%; Mouse: 100%; Bovine: 100%; Rabbit: 100%; Guinea pig: 100%; Zebrafish: 90%; Horse: 75%
Protein Families:	Transcription Factors



[View online »](#)

Product images:

WB Suggested Anti-ZNF449 Antibody Titration:
0.2-1 ug/ml; Positive Control: Transfected 293T